

Nitro Aluminum Rollator, 10" Casters



HCPCS:




FIGURE A



FIGURE B




FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G

- Attractive, Euro-style design
- Brake cable inside frame for added safety
- Handle height easily adjusts with unique push button (Figure A)
- Back support height easily adjusts with tool-free thumb screw (Figure B)
- Caster fork design enhances turning radius
- Large 10" front casters allow optimal steering and rolling comfort
- Easily folds with one hand to ultra-compact size for storage (Figure C)
- Handle is located on seat
- Plastic clip holds the Nitro closed when frame is folded (Figure D)
- Cross-brace design allows for side-to-side folding and added stability
- Seat is durable and comfortable (Figure E)
- Handsome and removable zippered storage bag Unique bag attachment keeps bag securely in place when rollator is open or folded (Figure F)
- Lightweight, aluminum frame

PRODUCTS

Description

RTL10266	Nitro Rollator, Red Frame, 1/ea
RTL10266BK	Nitro Rollator, Black frame,1EA
RTL10266WT	Nitro Rollator, White frame,1EA

SPECIFICATIONS

Description

5 Year Limited Warranty	(brake and handle housing)
Carton Shipping Weight	21.5 lbs.
Handle Height	33.5" - 38.25"
Seat Dimensions	18" (W) x 20.5" (H) x 10" (D)
Unit Dimensions	23" (W) x 27.75" (L)
Warranty Limited Lifetime	(frame)
Weight	17.5 lbs.
Weight Capacity	300 lbs

ACCESSORIES

Description

892	Cane Holder, Nitro Rollator, 1/ea, EA
10266-CH	Nitro Cup Holder, 1/ea, EA

PARTS

Description

1026603	HandGripRight for Nitro,1/ea, EA
1026602	HandGripLeft for Nitro,1/ea, EA

PARTS	Description
1026601-RN	HandBrakeAssembly,Left,10266,1ea, EA
1026601-G	HandBrakeAssembly,Left,10266BK/10266BK-T, EA
1026601-W	HandBrakeAssembly,Left,10266WT,1/ea, EA
1026630N	HandlebarForNitro/NitroTall,Side Mount, EA
1026604-B	BackrestAssemblyfor 10266 & 10266WT,1/ea, EA
1026605-W	BackrestLockingKnobW/Boltfor10266WT,1/pr, PR
1026605-G	BackrestLockingKnobW/Boltfor10266BK,1/pr, PR
1026604-G	Backrest Assembly for 10266BK,1/ea, EA
1026605-R	BackrestLockingKnobW/Boltfor10266,1/pr, PR
1026624	RightForkCompleteW/CasterforNitro, 1/ea, EA
1026625	LeftForkCompleteW/CasterforNitro, 1/ea, EA
1026617	FrontWheelWithBearingsforNitro, 1/ea, EA
1026616	WheelBrakeRightW/Hdwr,Nitro,Metal,1/ea, EA
1026615	WheelBrakeLeftW/Hdwr,Nitro,metal,1/ea, EA
1026612-R	RearWheelWithBearingsforNitro,Right,1/ea, EA
1026612-L	RearWheelWithBearingsforNitro,Left,1/ea, EA
1026609-FSW	SeatLiftingStrapW/AttachHrdwrefor10266WT, EA
1026609-FSG	SeatLiftingStrapW/AttachHrdwrefor10266BK, EA