

Nitro® Glide Knee Walker

HCPCS:



The Nitro Knee Walker offers users recovering from a below-the-knee injury the ability to remain mobile and glide through their daily activities in an all-new sleek design. Offered in two vibrant colors, this patented design features width expanding front axels to provide the user with adjustability and support to meet desired turning radius when navigating in different indoor and outdoor settings. The ergonomic design of the adjustable knee pad provides the necessary comfort and elevation to the injured leg while the adjustable handlebar also accommodates to user preference. The lightweight frame safely encloses the brake cable inside to enhance its modern design and the foldable feature makes it convenient for storage and transport.

- Key Features & Benefits
- Width expanding front axels (16", 18", 20") for adjustment of turning radius
- Easy-to-replace brake cables enclosed in frame tubing
- Tool-free height adjustment of handle (30"- 39") and knee pad (17" - 23").
- Offset, reversible knee pad for left or right leg use
- Includes removable storage basket with handle
- Lightweight and foldable frame for convenient storage and transport
- Simple-to-operate hand brakes control dual rear braking system

PRODUCTS	Description
791BL	Nitro Knee Walker, Blue, 1/EA
791RD	Nitro Knee Walker, Red, 1/EA
SPECIFICATIONS	Description
Warranty	Lifetime Limited Warranty
Product Weight	23 lb
Material	Aluminum, Plastic, Rubber, Foam
Handle Height	Min: 30" / Max: 39"
Weight Capacity	350 lb
Carton Shipping Weight	28 lb
Dimensions Assembled	35.3" (L) X 15.6" - 20.3" (W) X 30.7" - 38.6" (H)
Knee Pad Dimensions	13" (W) X 7" (D)
Knee Pad Height	Min: 17" / Max: 23"
CA Prop 65	no