

Phoenix LT 4-Wheel Scooter

HCPCS:




FIGURE A




FIGURE B



FIGURE C



FIGURE D



FIGURE E



FIGURE F



FIGURE G



FIGURE H

The Phoenix 4-wheel travel scooter combines a 350 lbs weight capacity with the take-apart portability of a travel scooter. The 4-wheel design and anti-tip wheels assist with overall stability for a more secure experience. This scooter also has a height-adjustable seat with flip-back and adjustable-width arms to accommodate varying body types. The seat also helps with portability by folding flat and allowing for easier transporting or storage. Your Phoenix 4-Wheel scooter can easily and quickly disassemble into five pieces.

- Interchangeable color panels
- Seat swivels 180 degrees for easy transfer
- Headlight for increased visibility
- Breaks down to 5 pieces for easy transport
- Comes with rugged 9"x3" wheels
- Padded Delta tiller

PRODUCTS	Description
PHOENIXLT4	Phoenix LT 4-Wheel Travel POV 18" Seat
PHOENIXLT4-20	Phoenix LT 4-Wheel Travel POV 20" Seat
SPECIFICATIONS	Description
Warranty	24 Month Limited Warranty on Electronics
Product Weight	127 lb
Product Dimensions	42.3" (L) x 20.1" (W) x 36.6" (H)
Weight Capacity	350 lb
Top Speed	4 mph
Batteries (pair)	12V x 22AH
Battery Charger	2A Offboard
Climbing Angle	9°
Color	Gray, Red, Blue
Controller	PG 45A
Ground Clearance	3.9"
HCPCS	K0801
Maximum Range	14 miles
Motor	24V x 250W x 170
Seat Type	Stadium
Turning Radius	70"
Warranty on Batteries	12 Month Warranty on Batteries
Warranty on Frame	Limited Lifetime Warranty on Frame

SPECIFICATIONS**Description**

Wheel Size

Front: 9" x 3" | Rear: 9" x 3"